

**Instrument: Gallios**

Laser	Detector	Filter	Mirror	Some common fluorochromes
Violet - 405 nm	FL9	450/50	480SP	BV421, Pacific Blue, Zombie Violet, DAPI
	FL10	550/40	—	BV510, Pacific Orange, Zombie Aqua
Blue - 488 nm	FL1	525/40	550SP	BB515, AF488, FITC, GFP, Zombie Green
	FL2	575/30	595SP	PE
	FL3	620/30	655SP	PE-Texas Red, PI
	FL4	670LP	730SP	BB700, PerCP, PerCP-Cy5.5, PE-Cy5, PE-Cy5.5, 7-AAD
	FL5	—	755LP	Pe-Cy7
Red - 640 nm	FL6	660/20	710SP	AF647, APC
	FL7	725/20	750SP	AF700, APC-AF700
	FL8	—	755LP	APC-Cy7, APC-H7, APC-AF750